

GK+ Blood Glucose & β -Ketone Meter

METER SPECIFICATIONS



Test Both Glucose and Ketone Levels



No Code Key Required!



Exceeds All Accuracy Standards



Bluetooth Integrated




Store 1000 Readings



FDA Approved



Lifetime Warranty



FREE Keto-Mojo App Available with Unlimited Storage!

Visit Keto-Mojo.com for more information!



INSERT STRIP

STRIP PORT

Insert a test strip in the strip port at the top of the meter. The strip port light turns on automatically when the meter is on.

SIDE BUTTONS

Toggle between various functions and settings.

DISPLAY

Shows test results and related information.

MAIN BUTTON

Press to turn meter on/off. Use to switch between functions and to set settings.

HYGIENIC EJECTOR

Located on the back of your meter, press up on the ejector to remove the test strip safely.



GLUCOSE

Measurement Range: 20 - 600 mg/dL

KETONE

Measurement Range: 0.1 - 8.0 mmol/L

